BIJSLUITER

Lacrisan®

kunstraan

SAMENSTELLING

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artificial tears

COMPOSITION

Lacrisan contains a cellulose derivative and polyvidone. Both substances increase the viscosity of the tear film.

APPLICATION

Lacrisan is used in dogs and cats for all instances in which too little tear fluid is produced or when there is a risk of dry eyes, for example with general anesthesia. The drops can also be used for eye examinations (gonioscopy), to clean the eyes and to prevent and remove tear stains.

METHOD OF USE AND ROUTE OF ADMINISTRATION

The drops are applied locally to the eye using the dropper bottle. It is important that the tip of the dropper does not come into contact with the eye or the hairs around the eye.

To keep the eyes moist, when tear production is insufficient: apply 1 or more drops to the eye up to 1 drop every hour.

Moistening the eyes when the eyes are not closed during general anesthesia: apply a few drops to the eye.

Cleaning the eyes: apply 1 or more drops to the eye.

Preventing tear stains: apply 1 or 2 drops to the eye 2 to 3 times a day.

To remove tear stains: drop Lacrisan on a cotton ball and clean the area around the eye.

STORAGE CONDITIONS

Use within 30 days after opening.

Store in a dark place at room temperature (15°C - 25°C).

Keep out of the sight and reach of children.

PACKAGE

Dropper bottle with 10 ml

MANUFACTURER

Produlab Pharma B.V.

Forellenweg 16

4941 SJ Raamsdonksveer

The Netherlands

DISTRIBUTOR

Alfasan Diergeneesmiddelen B.V.

Kuipersweg 9

3449 JA Woerden

The Netherlands